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for sterilizable silicone products

 Instructions for use

This reprocessing instructions can be used for the following fitting sets:

- MED1001770 Cramer Pessary Fitting Set
- MED1000489 Ring Pessary Fitting Set
- MED1001834 Ring Pessary thick Fitting Set
- MED1000656 Sieve Bowl Pessary Fitting Set
- MED1001835 Sieve Bowl Pessary EF Fitting Set
- MED1001733 T-Pessary Fitting Set Type B
- MED1001732 T-Pessary Fitting Set Type R
- MED1000765 Urethra Pessary Fitting Set
- MED1000488 Urethra Bowl Pessary Fitting Set
- MED1000487 Cube Pessary Fitting Set
  
- MED1001868 Vaginal Phantom 30 mm × 240 mm
- MED1001833 Vaginal Phantom 32 mm × (2 × 150 mm)
  
- MED1000475 Diaphragm Fitting Set
  
- MED1000367 Mouth Wedge for adults, 20 × 25 × 90 mm, blue, standard
- MED1000368 Mouth Wedge for adults, 20 × 25 × 90 mm, transparent, soft
- MED1000658 Mouth Wedge for children, 17 × 14 × 80 mm, blue, standard
- MED1000657 Mouth Wedge for children, 17 × 14 × 80 mm, transparent, soft
- MED1001086 Mouth Wedge wedge-shaped, 90 × 35 mm, blue, standard
- MED1001087 Mouth Wedge wedge-shaped, 90 × 35 mm, transparent, soft

## Important general safety information:

medesign products are generally delivered in a non-sterile state and must be cleaned and sterilised according to these instructions prior to first use.

Please read these instructions for reprocessing silicone products thoroughly, entirely and carefully before using the product for the first time. Always follow these instructions and save them for further reference.

Observe the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM) on 'hygiene requirements in the reprocessing of medical devices'.

This product may not be reused and has to be destroyed if it has been used on a patient with Creutzfeldt-Jakob Disease (CJD) or someone at risk of CJD.

Only professionally qualified personnel may clean and sterilise silicone products.

In case of damage the silicone product has to be replaced immediately, also during use. Further use can lead to complications and / or endanger the user.

Hot-air sterilisation is prohibited as it destroys the silicone product.

Never clean silicone products with metal brushes or steel wool.

All products need to be removed from their transport packaging prior to use. Other packaging, protective covers, foil etc. has to be removed and disposed of in an environmentally responsible manner.

Material durability: Please ensure that the cleaning and disinfecting agents do not contain any of the following ingredients: strong acids (pH < 5 / oxidising acids), alkalis (> pH 10), organic solvents, alcohol, phenol, ammonia, petrol, halogens, alkyl halides, higher concentrations of sodium chloride, oxidising agents.

Maintenance and repair may solely be carried out by the company medesign I.C. GmbH or by persons authorised by medesign I.C. GmbH.

Violation leads to invalidity of guarantee and warranty claims.

In order to protect patients, users and third parties all products have to be cleaned, disinfected and sterilised according to this instruction manual prior to each utilisation and before returning goods. This specifically also applies to the first use after delivery, as all products and packaging are delivered non-sterile. Effective cleaning and disinfection is an essential requirement for successful sterilisation. Please ensure, within the scope of your responsibility for the sterility of the products in use, that only sufficient device- and product-specific validated methods are used for cleaning / disinfection and sterilisation. Please have the used devices, such as the disinfectant or steriliser, serviced and checked regularly and ensure that the validated parameters are adhered to with every cycle.

In addition observe your country's specific rules and regulations concerning medical practices or hospitals.

**1. Indications and field of application:**

To be used by adequately trained and qualified healthcare professionals in various medical fields. Further information is to be found in the instruction manuals of the individual silicone products.

**2. Contraindications:**

CJK – Creutzfeldt-Jakob disease

BSE – Bovine spongiform encephalopathy; so-called „mad cow disease“ (e.g. Creutzfeldt-Jakob disease)

TSE – Transmissible spongiform encephalopathy

Based on the patient's general state of health the responsible physician has to decide whether the planned procedure can be done. Further indications are to be found in relevant specialist literature. In the case of a suspected diagnosis of Creutzfeldt-Jakob disease measures to prevent a possible transmission to other patients, users or third parties have to be taken. Country-specific regulations concerning preparation have to be strictly adhered to.

**3. Application:**

According to the individual instruction manuals.

**4. Inspection:**

Important!

Inspections are to be done prior to and after each utilisation. Discontinue use of silicone products that are damaged, incomplete or have loose parts. Do not attempt to repair on your own.

Check products for damage, loose or missing parts and rough surfaces. Hot-air sterilisation is prohibited as it destroys silicone.

**In addition please adhere to the points in 5.7 Inspection/Functional check/Care.**

## **5. Reprocessing procedures as described in chapter 5.1 and 5.11: alkaline, manual ultrasonic precleaning, mechanical cleaning and disinfection using a device**

### **5.1 Limitations and restrictions with regard to reprocessing:**

Due to the product design and the materials used it is not possible to determine a defined limit regarding the maximum number of reprocessing cycles. The end of the product's life cycle is usually established based on wear and tear and damage during use.

### **5.2 Preparation at place of use:**

The silicone products should preferably be cleaned and disinfected immediately after use.

Remove visible impurities from the silicone product directly after use. These impurities should not dry out on the object, as this will complicate disinfection and cleaning further.

Do not use any fixing agents or hot water ( $< 40\text{ }^{\circ}\text{C}$ ), as this will lead to fixing of residues, which in turn can influence successful cleaning.

### **5.3 Transportation:**

Safe storage of the silicone products in a sealed container during transportation to the place of preparation is paramount to avoid damage to and contamination of the environment.

### **5.4 Cleaning:**

Place silicone products in cold water for at least 5 minutes. If possible clean the silicone products under running cold water with a soft brush until no more residue is visible.

In the case of cavities, hollows or drilled openings pressure spray with a water spray for at least 10 seconds (pulse technique).

Where possible place the silicone product in an ultrasonic bath containing demineralised water at  $40\text{ }^{\circ}\text{C}$  with a 0.50 % alkaline cleaning agent and treat with ultrasound. The ultrasound frequency should be at least 35 kHz. Remove instruments and rinse in lukewarm water ( $< 40\text{ }^{\circ}\text{C}$ ).

The cleaning performance must be checked by visual inspection. The silicone products must be optically spotlessly clean, if this is not the case the process must be repeated. The inspection is done visually, critical areas like cavities, hollows and drilled openings require special, thorough checking.

Fresh disinfecting and cleaning agents must be used daily. By using cleaning agents for an extended length of time the following problems can arise:

- Reduced effect regarding disinfection due to pollution by impurities.

The residues resulting from the cleaning process must definitely be removed otherwise stains and discolouration will appear on the silicone products.

## **5.5 Disinfection:**

The national compliances regarding the A0 values (see ISO 15883) must be complied with when cleaning using a machine with thermal disinfection.

Please note:

- The sieve tray may not be overloaded so that all instruments can be thoroughly rinsed. Instruments with long or narrow cavities or hollows must also be perfused internally.

## **5.6 Drying:**

The exterior of the silicone product is dried during the drying cycle of the cleaning / disinfecting appliance. If necessary, further manual drying, using a lint-free cloth, is possible. Hollows and cavities of silicone products must be cleaned with sterile compressed air.

## **5.7 Inspection / Functional check / Care:**

Visually assess the cleanliness of the silicone products. If necessary repeat reprocessing until the silicone product is optically clean. Care as described below and functional check as described in paragraph 5 to be implemented.

Stained silicone products are the result of defects in the reprocessing. Causes of staining could possibly be due to:

- inadequate mechanical or manual cleaning
- unsuitable cleaning agents, disinfection agents and silicone care products
- non-compliance with dosage specifications of the cleaning agents, disinfection agents and silicone care products or residues (carry over)
- water influences e.g. unfamiliar ions such as iron and silicate
- residues from pharmaceuticals, markers or chemical indicators
- procedural errors

Check all products for damaged surfaces, splitting, distortion, flexibility / function, cracking etc. after cleaning i.e. cleaning and disinfection and discard damaged or otherwise worn silicone products.

Please forward only reprocessed products (cleaned, disinfected and sterilised) to the manufacturer for repair.

**5.8 Packaging:**

We recommend that sterilisation be carried out in disposable sterilisation packaging and / or sterilisation containers. These must comply with the legal requirements (e.g. DIN EN ISO / ANSI AAMI ISO 11607 / ISO 11607 / EN 868), must be suitable for the sterilisation process and offer sufficient protection of products i.e. prevent mechanical damage to sterilisation packaging.

Please perform regular maintenance services as recommended by the manufacturer.

**5.9 Sterilisation:**

Sterilisation of the products using a fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) taking into account the respective national requirements:

- pre-vacuum phases with a minimum pressure of 60 millibars
- heating at a sterilisation temperature of at least 132 °C; max. 137 °C
- shortest holding time: 4 minutes
- drying time: at least 10 minutes

Please take note of the following:

The sterilisation equipment must comply with the applicable CEN- i.e. DIN-norms.

The instruction manual must be followed unconditionally.

**5.10 Storage:**

All sterilised products must be stored in a clean, dry, dust-free place away from direct sunlight and at a moderate temperature between 5 and 40 °C. Transportation and storage may not adversely affect the properties of the medical device.

- relative humidity 30 – 70 %
- protect from intensive sunlight and artificial light

**5.11 Information for validating the reprocessing:**

The following testing instructions, materials and machines were used for the validation:

Cleaning agent:	Neodisher MediClean forte; Dr. Weigert; Hamburg
Neutraliser:	Neodisher Z; Dr. Weigert; Hamburg
Cleaning and disinfection device:	Miele
Slide-in module:	Miele injector unit E450 / 1

Additional instructions: Should the chemicals / machines described above not be available it is the responsibility of the user to validate the process used. It is obligatory for the user to ensure that the reprocessing procedure, including resources, material and personnel are suitable to achieve the results required. The best available technology and national laws demand the compliance with validated processes.

**Disposal:**

Used or damaged silicone products can be disposed of when in a sterile condition. The country-specific legal requirements and regulations must be complied with when disposing of silicone products.

**Warranty:**

We warrant a full guarantee in the case of quality or production defects. In the event of obvious shortcomings, which have arisen due to manufacturing errors or the use of defective materials the products will be processed or replaced at no extra cost. The claim under guarantee is excluded in the case of damage due to inappropriate handling e.g. mechanical influence, dropping, overload etc. The claim under guarantee expires in the case of repairs by unauthorised people.

**Applied standards and regulations:**

This product conforms with Regulation 93 / 42 / EWG Class I, IIa und IIb. Details may be found in the respective instruction manuals.

Manufacturers symbol

Symbol for unsterile

Follow directions for use

Symbol for store in a dry place

Identification of products introduced in accordance with regulation 93 / 42 EWG

Protect from direct sunlight

Batch description / batch number

Please note: Attention should be paid to enclosed documents, observe the instruction manual

## Legal requirements, regulations and standards:

93 / 42 / EWG	Medical Devices Directive
MPG	Medical Device Act
DIN EN 285	Sterilisation – steam sterilisers – large sterilisers
DIN EN 868	Packaging materials and systems for medical devices to be sterilised
DIN EN 13060	Small steam sterilisers
DIN EN ISO 11607	Packaging for medical devices to be sterilised in the final packaging
DIN EN ISO 14937	Sterilisation of products for medical care
DIN EN ISO 14971	Medical devices – application of risk management
DIN EN ISO 15883	Cleaning / Disinfection appliances (validation and operation)
DIN EN ISO 17664	Sterilisation of medical devices – Information to be provided by the manufacturer
DIN EN ISO 17665-1	Sterilisation of products for medical care

medesign I.C. GmbH  
Dietramszeller Straße 6  
83623 Dietramszell - Linden  
Germany

Tel.: +49 (0)8027 90738-0  
Fax: +49 (0)8027 90738-27  
www.medesign.de  
info@medesign.de

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